

Streptococcus Group A antibody (FITC)

Cat. No. GTX36330

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	ICC/IF, ELISA, Sandwich ELISA
Reactivity	Streptococcus Group A

Package
500 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Acetone fixation of the antigen source is recommended prior to staining.

Functions both as capture and reporter antibody.

Capture: GTX36558, Detection:GTX36330/GTX36482/GTX36406

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Streptococci, group A
Purification	Purified by antigen-affinity chromatography
Conjugation	Fluorescein isothiocyanate (FITC)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).