

# Cytomegalovirus gB antibody (HRP)

**Cat. No. GTX36525**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Application</b>	WB, ICC/IF, ELISA, IHC
<b>Reactivity</b>	Cytomegalovirus

**Package**  
500 µl

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:20-1:200
ICC/IF	Assay dependent
ELISA	1:200-1:1000
IHC	1:20-1:200

Not tested in other applications.

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 1% BSA
<b>Preservative</b>	0.002% Thimerosal
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1-2 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	purified virions of strain AD169
<b>Purification</b>	IgG fraction
<b>Conjugation</b>	Horseradish peroxidase(HRP)

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).