

## EBV gp 250/350 antibody [0221]

## Cat. No. GTX36648

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, ICC/IF
Reactivity	Epstein-Barr virus

Package 500 μg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:10-1:50

Not tested in other applications.

Store at 4°C.
gnostic use in humans or animals. Not for animal or
py, reverse engineer or otherwise attempt to
9



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 10 Page 1 of 1