

## Giardia lamblia 65 kDa antibody [6261]

## Cat. No. GTX36744

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2b
Applications	WB, ICC/IF, ELISA, Sandwich ELISA
Reactivity	Giardia lamblia

Package 500 μΙ

## **Applications**

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:10-1:50
ELISA	1:20-1:200
Sandwich ELISA	Assay dependent

Note: Capture: GTX03385, Detection: GTX36744

The detection antibody requires HRP conjugation. Please contact us for the conjugation service.

Not tested in other applications.

Properties		
Form	Liquid	
Buffer	PBS	
Preservative	0.1% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.	
Concentration	$100  \mu g/ml$ (Please refer to the vial label for the specific concentration.)	
Purification	Protein A purified	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 04 Page 1 of 1