

SARS-CoV Nucleoprotein antibody [3851]

Cat No. GTX36802

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	WB, ICC/IF, ELISA
Reactivity	SARS Coronavirus

Reference (1)
Package
500 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Dilution
WB	Assay dependent
ICC/IF	1:10-1:50
ELISA	1:20-1:200

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	0.01M PBS pH7.2
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).