

# West Nile virus Envelope antibody [3671]

## Cat. No. GTX36805

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2b
Applications	ICC/IF, ELISA
Reactivity	West Nile virus

Package 500 μg

## Applications

#### **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200

Not tested in other applications.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 $\mu$ g/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 06 Page 1 of 2



### DATA IMAGES



#### GTX36805 Image

West Nile Virus monoclonal antibody cross reactivities



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 06 Page 2 of 2