

## West Nile virus Envelope antibody [3655]

Cat. No. GTX36807

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, ELISA
Reactivity	West Nile virus

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200

Not tested in other applications.

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

Virus Reactivity (IFA)											
	*	*	*	*	*	*	*	*	*	*	*
Maby	WNV	SLE	JEV	YFV	DEN	TBE	EEE	WEE	POW	LAC	VEE
3631	pos	pos	neg	neg	neg	neg	neg	neg	neg	neg	neg
3641	pos	pos	neg	neg	neg	neg	neg	neg	neg	neg	neg
3655	pos	pos	neg	pos	neg	neg	neg	neg	neg	neg	neg
3661	pos	pos	neg	pos	neg	neg	neg	neg	neg	neg	neg
3671	pos	neg	pos	neg	neg	pos	neg	neg	neg	neg	neg

\*Flaviviruses

**GTX36807 Image**

West Nile Virus monoclonal antibody cross reactivities



For full product information, images and publications, please visit our [website](#).