

Hepatitis B virus Surface Antigen antibody [1837]

Cat. No. GTX36858

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, IHC, Sandwich ELISA
Reactivity	Hepatitis B virus

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:20-1:200
IHC	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX36859/GTX36860/GTX36861, Detection : GTX36858.
The detection antibody requires HRP conjugation. Please contact us for the conjugation service.**

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native antigen
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).