

Hepatitis B virus Surface Antigen antibody [1836]

Cat. No. GTX36859

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, IHC, Sandwich ELISA
Reactivity	Hepatitis B virus

References (1)

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:20-1:200
IHC	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX36859, Detection : GTX36858/GTX40878.

Please notice that GTX36858 need to be conjugated to HRP to function as the detection antibody when paired with GTX36859. Please contact us for conjugation service.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native antigen
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 04 / 06 Page 1 of 1