

## Hepatitis B virus Surface Antigen antibody [1835]

**Cat. No. GTX36860**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	WB, ELISA, IHC, Sandwich ELISA
<b>Reactivity</b>	Hepatitis B virus

Package  
500 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:20-1:200
IHC	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX36860, Detection : GTX36858/GTX40878.**

**Please notice that GTX36858 need to be conjugated to HRP to function as the detection antibody when paired with GTX36860. Please contact us for conjugation service.**

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	100 µg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Native antigen
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).