

Chlamydia trachomatis LPS antibody [1681]

Cat. No. GTX36874

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, ELISA, IHC
Reactivity	Chlamydia trachomatis, Chlamydia pneumoniae, Chlamydia psittacii

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200
IHC	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX36874 ICC/IF Image

ICC/IF analysis of L929 cells with Chlamydia trachomatis LGV Type 1 infection using GTX36874 Chlamydia trachomatis LPS antibody [1681].

Red : Uninfected cells

Green : Primary antibody; Chlamydia elementary bodies and inclusions of the same.



For full product information, images and publications, please visit our [website](#).