

# alpha 1 Antichymotrypsin antibody [801]

**Cat. No. GTX36929**

|                   |            |
|-------------------|------------|
| <b>Host</b>       | Mouse      |
| <b>Clonality</b>  | Monoclonal |
| <b>Isotype</b>    | IgG2a      |
| <b>Reactivity</b> | Human      |

**Package**  
100 µg

## Properties

|                      |                                                                                                                                                                                                                                                                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                                                                                                                                         |
| <b>Preservative</b>  | No preservatives                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| <b>Concentration</b> | 1.30 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                 |
| <b>Immunogen</b>     | Full length protein: Human Serum Alpha-1-Antichymotrypsin                                                                                                                                                                                                                                                                                   |
| <b>Purification</b>  | Protein G purified                                                                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| <b>Note</b>          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).