

GSK3 beta antibody

Cat. No. GTX37303

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100-1:200
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized non-phosphopeptide derived from human GSK3 beta around the phosphorylation site of serine 9 (T -T-S(p)-F-A).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX37303 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX37303 GSK3 beta antibody.



GTX37303 WB Image

WB analysis of A2780 and HT-29 cell lysates using GTX37303 GSK3 beta antibody.



GTX37303 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX37303 GSK3 beta antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).