

MCP1 / CCL2 antibody

Cat. No. GTX37379

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rabbit

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:500

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 34-89 amino acids of human MCP1
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX37379 IHC-P Image


IHC-P analysis of human rectal carcinoma tissue using GTX37379 MCP1 / CCL2 antibody.

Dilution : 1:200


GTX37379 IHC-P Image

IHC-P analysis of mouse bone tissue using GTX37379 MCP1 / CCL2 antibody.

Dilution : 1:150


GTX37379 IHC-P Image

IHC-P analysis of rabbit aorta tissue using GTX37379 MCP1 / CCL2 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](#).