

LEF1 antibody

Cat. No. GTX37387

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, FCM
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 341-389 amino acids of human LEF1
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX37387 FCM Image

FCM analysis of 4% PFA fixed Molt4 cells using GTX37387 LEF1 antibody.


Green : Primary antibody

Orange : Isotype control

Blue : Negative control

Permeabilization : 90% ice-cold methanol


Dilution : 1 µg



GTX37387 WB Image

WB analysis of various mouse tissue lysates using GTX37387 LEF1 antibody.

Dilution : 1:1000



GTX37387 WB Image

WB analysis of mouse lymph node tissue lysate using GTX37387 LEF1 antibody.

Dilution : 1:300

Loading : 30 µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).