

IL11 antibody

Cat. No. GTX37407

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100-2000
IHC-P	1:100-800 (based on 0.5 mg/ml)

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human IL11.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX37407 WB Image

WB analysis of rat liver (Lane 1), mouse heart (Lane 2), rat colon (Lane 3), mouse spleen (Lane 4), mouse kidney (Lane 5), rat testis (Lane 6) lysates using GTX37407 IL11 antibody.


Dilution : 1µg/ml



GTX37407 IHC-P Image

IHC-P analysis of rat colon tissue using GTX37407 IL11 antibody.


Dilution : 2.5µg/ml



GTX37407 IHC-P Image

IHC-P analysis of rat skin tissue using GTX37407 IL11 antibody.

Dilution : 2.5µg/ml



GTX37407 IHC-P Image

IHC-P analysis of rat skin tissue using GTX37407 IL11 antibody.

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).