

FOXJ1 antibody

Cat. No. GTX37424

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
IHC-P	1:400-800

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 171-250 amino acids of human FOXJ1
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX37424 WB Image


WB analysis of 40µg of mouse and rat testis tissue lysates using GTX37424 FOXJ1 antibody.


GTX37424 IHC-P Image

IHC-P analysis of human brain glioma tissue using GTX37424 FOXJ1 antibody.


GTX37424 IHC-P Image

IHC-P analysis of rat brain tissue using GTX37424 FOXJ1 antibody.


GTX37424 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX37424 FOXJ1 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).