

Endothelin 1 antibody

Cat. No. GTX37427

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100-1000
ICC/IF	1:200

Not tested in other applications.

Calculated MW 24 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Endothelin 1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX37427 ICC/IF Image


ICC/IF analysis of HeLa cells using GTX37427 Endothelin 1 antibody.

Dilution : 2.5µg/ml


GTX37427 WB Image

WB analysis of HeLa cells (Lane 1), mouse lung (Lane 2), rat heart (Lane 3), rat large intestines (Lane 4) lysates using GTX37427 Endothelin 1 antibody.

Dilution : 1µg/ml


GTX37427 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX37427 Endothelin 1 antibody.

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).