

Elastin antibody

Cat. No. GTX37428

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat

References (4)

★★★★★ Review (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:100-500
IHC-Fr	Assay dependent

Not tested in other applications.

Calculated MW 68 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Elastin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX37428 WB Image

WB analysis of mouse spinal cord (Lane 1), rat bladder (Lane 2), mouse stomach (Lane 3), rat lung (Lane 4), rat testis (Lane 5) tissue lysates using GTX37428 Elastin antibody.


Dilution : 1µg/ml



GTX37428 IHC-P Image

IHC-P analysis of rat spinal cord tissue using GTX37428 Elastin antibody.


Dilution : 2.5µg/ml



GTX37428 IHC-P Image

IHC-P analysis of mouse stomach tissue using GTX37428 Elastin antibody.


Dilution : 2.5µg/ml



GTX37428 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX37428 Elastin antibody.

Dilution : 1:100



GTX37428 IHC-Fr Image

IHC-Fr analysis of mouse skin tissue using GTX37428 Elastin antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).