

EAAT2 antibody

Cat. No. GTX37432

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 111-190 amino acids of human EAAT2
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX37432 IHC-P Image

IHC-P analysis of rat brain tissue using GTX37432 EAAT2 antibody.

Fixation : 4% PFA


Antigen retrieval : boiling with citrate buffer

Dilution : 1:200


GTX37432 IHC-P Image

IHC-P analysis of rat brain tissue using GTX37432 EAAT2 antibody.

Dilution : 1:200


GTX37432 IHC-P Image

IHC-P analysis of rat brain tissue using GTX37432 EAAT2 antibody.

Fixation : 4% PFA

Antigen retrieval : boiling with citrate buffer

Dilution : 1:200



For full product information, images and publications, please visit our [website](#).