

# NKp44 antibody

**Cat. No. GTX37443**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human, Rat |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| IHC-P              | Assay dependent      |

Not tested in other applications.

**Calculated MW** 31 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | 10mM TBS, 0.5% BSA, 25% Glycerol                                                                                                                                                                                           |
| <b>Preservative</b>  | 0.015% ProClin 300                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| <b>Immunogen</b>     | KLH conjugated synthetic peptide derived between 211-266 amino acids of human NKp44                                                                                                                                        |
| <b>Purification</b>  | Protein A purified                                                                                                                                                                                                         |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


## DATA IMAGES


**GTX37443 IHC-P Image**

IHC-P analysis of human laryngeal cancer tissue using GTX37443 NKp44 antibody.


**GTX37443 IHC-P Image**

IHC-P analysis of human laryngeal cancer tissue using GTX37443 NKp44 antibody.


**GTX37443 WB Image**

WB analysis of rat liver tissue lysate using GTX37443 NKp44 antibody.



For full product information, images and publications, please visit our [website](#).