

CD5L antibody

Cat. No. GTX37448

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM, IHC
Reactivity	Mouse, Rat

References (3)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	1:400-800
FCM	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	0.01M TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 250-300 amino acids of human CD5L
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx37448 IHC-P Image


IHC-P analysis of rat spleen tissue using GTx37448 CD5L antibody.
Antigen retrieval : Boiling in sodium citrate buffer (pH6.0) for 15min
Dilution : 1:400


GTx37448 IHC-P Image

IHC-P analysis of mouse lymph tissue using GTx37448 CD5L antibody.


GTx37448 IHC Image

The data was published in the journal Nat Commun in 2018. [PMID: 29500407](https://pubmed.ncbi.nlm.nih.gov/29500407/)


GTx37448 ICC/IF Image

The data was published in the journal Nat Commun in 2018. [PMID: 29500407](https://pubmed.ncbi.nlm.nih.gov/29500407/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).