

CD35 antibody

Cat. No. GTX37450

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100-2000
IHC-P	1:100-800 (based on 0.5 mg/ml)

Not tested in other applications.

Calculated MW 224 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human CD35.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX37450 WB Image


WB analysis of rat pancreas (Lane 1), rat thymus (Lane 2), human breast cancer (Lane 3) tissue lysates using GTX37450 CD35 antibody.

Dilution : 1µg/ml


GTX37450 IHC-P Image


IHC-P analysis of mouse lymph node tissue using GTX37450 CD35 antibody.

Dilution : 2.5µg/ml


GTX37450 IHC-P Image

IHC-P analysis of human endometrial cancer tissue using GTX37450 CD35 antibody.

Dilution : 2.5µg/ml


GTX37450 IHC-P Image

IHC-P analysis of rat stomach tissue using GTX37450 CD35 antibody.

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).