

Reelin antibody

Cat. No. GTX37552

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-800 (based on 0.5 mg/ml)
IHC-Fr	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Reelin. Please contact us for the exact immunogen sequence.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX37552 IHC-P Image

IHC-P analysis of human mammary fibroma tissue using GTX37552 Reelin antibody.
Dilution : 2.5µg/ml




GTX37552 IHC-P Image

IHC-P analysis of mouse skin tissue using GTX37552 Reelin antibody.
Dilution : 2.5µg/ml



GTX37552 IHC-Fr Image

IHC-Fr analysis of mouse brain tissue using GTX37552 Reelin antibody.
Dilution : 2.5µg/ml



GTX37552 IHC-P Image

IHC-P analysis of human endometrial cancer tissue using GTX37552 Reelin antibody.
Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](#).