

Nociceptin Receptor antibody

Cat. No. GTX37594

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
IHC-P	1:400-800
IHC-Fr	Assay dependent

Not tested in other applications.

Calculated MW 41 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived between 13820 amino acids of human Nociceptin Receptor
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx37594 IHC-P Image


IHC-P analysis of rat brain tissue using GTx37594 Nociceptin Receptor antibody.


GTx37594 IHC-Fr Image

IHC-Fr analysis of human placenta tissue using GTx37594 Nociceptin Receptor antibody.


GTx37594 WB Image

WB analysis of rat heart (Lane 1) and rat brain(Lane 2) tissue lysates using GTx37594 Nociceptin Receptor antibody.


GTx37594 IHC-Fr Image

IHC-Fr analysis of mouse hippocampus tissue using GTx37594 Nociceptin Receptor antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).