

## Fibronectin antibody

**Cat. No. GTX37694**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Reactivity</b>	Mouse, Rat

References ( 1 )

Package

100 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:400-800
IHC-Fr	Assay dependent

Not tested in other applications.

**Calculated MW** 263 kDa. ([Note](#))

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	10mM TBS, 0.5% BSA, 25% Glycerol
<b>Preservative</b>	0.015% ProClin 300
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	KLH conjugated synthetic peptide derived between 1211-1290 amino acids of human Fibronectin
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated


**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).


## DATA IMAGES

**GTX37694 IHC-Fr Image**

IHC-Fr analysis of rat pancreas tissue using GTX37694 Fibronectin antibody.

**GTX37694 IHC-P Image**

IHC-P analysis of rat kidney tissue using GTX37694 Fibronectin antibody.

**GTX37694 IHC-P Image**

IHC-P analysis of rat kidney tissue using GTX37694 Fibronectin antibody.



For full product information, images and publications, please visit our [website](#).