

beta 2 Defensin antibody

Cat. No. GTX37717

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100-1500
IHC-P	1:100-500

Not tested in other applications.

Calculated MW 8 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human beta 2 Defensin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX37717 IHC-P Image

IHC-P analysis of rat trachea tissue using GTX37717 beta 2 Defensin antibody.
Dilution : 2.5µg/ml


GTX37717 IHC-P Image

IHC-P analysis of rat lung tissue using GTX37717 beta 2 Defensin antibody.
Dilution : 2.5µg/ml


GTX37717 IHC-P Image

IHC-P analysis of rat uterus tissue using GTX37717 beta 2 Defensin antibody.
Dilution : 2.5µg/ml


GTX37717 WB Image

WB analysis of rat skin (Lane 1), rat trachea (Lane 2), mouse lung (Lane 3), rat kidney (Lane 4), mouse uterus (Lane 5), U251 (Lane 6) lysates using GTX37717 beta 2 Defensin antibody.
Dilution : 1µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).