

Collagen X antibody

Cat. No. GTX37732

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC
Reactivity	Human, Mouse, Rat, Rabbit

References (22)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:200-2000
IHC-P	1:300-500
IHC	Assay dependent

Not tested in other applications.

Calculated MW 66 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Collagen X.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx37732 IHC-P Image

IHC-P analysis of rat liver tissue using GTx37732 Collagen X antibody.


Dilution : 2.5µg/ml



GTx37732 WB Image

WB analysis of mouse pancreas (Lane 1), human breast cancer (Lane 2), human mammary fibroma (Lane 3) lysates using GTx37732 Collagen X antibody.


Dilution : 1µg/ml



GTx37732 WB Image

WB analysis of mouse spinal cord (Lane 1), mouse liver (Lane 2), mouse pancreas (Lane 3), human breast cancer (Lane 4), human mammary fibroma (Lane 5) lysates using GTx37732 Collagen X antibody.

Dilution : 1µg/ml



GTx37732 IHC-P Image

IHC-P analysis of mouse bone marrow tissue using GTx37732 Collagen X antibody.

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).