

CD24 antibody

Cat. No. GTX37755

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Rat

References (11)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 8 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 0.5% BSA, 25% Glycerol
Preservative	0.015% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human CD24 Propeptide:60-80/80.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX37755 IHC-P Image

IHC-P analysis of rat lung tissue using GTX37755 CD24 antibody.
Dilution : 1:100




GTX37755 WB Image

WB analysis of human colon carcinoma (Lane 1) and rat heart (Lane 2) tissue lysates using GTX37755 CD24 antibody.



GTX37755 IHC-P Image

IHC-P analysis of human colon carcinoma tissue using GTX37755 CD24 antibody.



GTX37755 IHC-P Image

The data was published in the 2020 in Int J Mol Sci. [PMID: 33375383](https://pubmed.ncbi.nlm.nih.gov/33375383/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).