

Her3 / ErbB3 (phospho Tyr1328) antibody

Cat. No. GTX38571

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100-1:200
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW	148 kDa. (Note)
Product Note	This antibody detects endogenous levels of Her3/ErbB3 only when phosphorylated at tyrosine 1328.

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human Her3/ErbB3 around the phosphorylation site of tyrosine 1328 (P -D-Y(p)-W-H).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 09 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX38571 WB Image

WB analysis of A431 cell lysate using GTX38571 Her3 / ErbB3 (phospho Tyr1328) antibody.



GTX38571 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX38571 Her3 / ErbB3 (phospho Tyr1328) antibody.



GTX38571 ICC/IF Image

ICC/IF analysis of MCF-7 cells using GTX38571 Her3 / ErbB3 (phospho Tyr1328) antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).