

p38 MAPK (phospho Thr180) antibody

Cat. No. GTX38595

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 41 kDa. ([Note](#))

Product Note p38 MAPK (phospho-Thr180) antibody detects endogenous levels of p38 MAPK only when phosphorylated at threonine 180.

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human p38 MAPK around the phosphorylation site of threonine 180 (E -M-T(p)-G-Y).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX38595 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX38595 p38 MAPK (phospho Thr180) antibody.


GTX38595 WB Image

WB analysis of HT-29 cell lysate using GTX38595 p38 MAPK (phospho Thr180) antibody.



For full product information, images and publications, please visit our [website](#).