

Stathmin 1 (phospho Ser25) antibody

Cat. No. GTX38657

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 17 kDa. ([Note](#))

Product Note This antibody detects endogenous levels of Stathmin 1 only when phosphorylated at serine 25.

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human Stathmin 1 around the phosphorylation site of serine 25 (I -L-S(p)-P-R).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX38657 WB Image

WB analysis of Jurkat cell lysate using GTX38657 Stathmin 1 (phospho Ser25) antibody.



GTX38657 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX38657 Stathmin 1 (phospho Ser25) antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).