

Vaccinia virus antibody [B608M]

Cat. No. GTX38723

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ICC/IF, ELISA
Reactivity	Vaccinia virus

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note Recognizes free Vaccinia virus and infected cells. Reactive with Lister and NYCBH strains.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native, intact Lister strain of Vaccinia virus.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).