

TSH beta antibody

Cat. No. GTX38745

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note < 0.1% cross reactivity with alpha hCG, LH, FSH, growth hormone, and prolactin.

Properties

Form	Liquid
Buffer	15mM Potassium Phosphate, 0.85% NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	4.90 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Intact TSH
Purification	Purified by Affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).