

Trypsin antibody [403]

Cat. No. GTX38777

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA, Immunoassay
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Immunoassay	Assay dependent

Not tested in other applications.

Product Note Does not cross-react with pancreatic chymotrypsin.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.90 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human pancreatic Trypsin.
Purification	Protein G purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).