

## V5 tag antibody [SV5-Pk3]

Cat. No. GTX38824

|              |                     |
|--------------|---------------------|
| Host         | Mouse               |
| Clonality    | Monoclonal          |
| Isotype      | IgG1                |
| Applications | WB, IP, ELISA       |
| Reactivity   | Species independent |

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1/100-1/500          |
| IP                 | Assay dependent      |
| ELISA              | Assay dependent      |

Not tested in other applications.

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | A recombinant peptide corresponding to the V5 tag fused to a class II - restricted T helper cell epitope.                                                                                                                  |
| Purification  | Protein G purified<br>From tissue culture supernatant                                                                                                                                                                      |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX38824 WB Image

WB analysis of V5 tagged protein using GTX38824 V5 tag antibody [SV5-Pk3].



For full product information, images and publications, please visit our [website](#).