

## Trypanosoma cruzi Lpg antibody [Tcr7D3]

Cat. No. GTX38864

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Application	ICC/IF, ELISA, Immunoblotting
Reactivity	Trypanosoma cruzi

Reference ( 1 )  
Package  
100 µl

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Immunoblotting	Assay dependent

Not tested in other applications.

## PROPERTIES

Form	Liquid
Buffer	Ascites
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	T. cruzi LPG purified from epimastigotes.
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).