

## Thyroxine (T4) antibody [057-11007]

**Cat. No. GTX38976**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	WB, EIA, Lateral Flow
<b>Reactivity</b>	Species independent

**Package**  
2 mg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
EIA	Assay dependent
Lateral Flow	Assay dependent

Not tested in other applications.

**Product Note** Cross-Reactivity:T3 <1% T2 <1% rT3 <1%

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	5.60 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	T4- BSA
<b>Purification</b>	Purified by Column Chromatography
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).