

Thyroxine (T4) antibody [057-11007]

Cat. No. GTX38976

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b

Package
2 mg

APPLICATION

Product Note Cross-Reactivity:T3 <1% T2 <1% rT3 <1%

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	5.60 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	T4- BSA
Purification	Purified by Column Chromatography
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).