

T Lymphocytes antibody [KEN-5]

Cat. No. GTX38978

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM
Reactivity	Rabbit

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent

Not tested in other applications.

Product Note

GTX38978 reacts with T cells. Recognises >90% thymocytes and 40-45% mesenteric lymph node cells. The antibody immunoprecipitates a 67kDa protein from thymocytes.

This antibody was initially reported as recognising rabbit CD5 (1), but does not bind to rabbit CD5 transfectants. Known rabbit CD5 antibodies also show binding to most B lymphocytes (2), which are not labelled by this clone.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Rabbit Thymocytes
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX38978 FCM Image

Staining of rabbit peripheral blood lymphocytes with Mouse anti-Rabbit T Lymphocytes (GTX38978)



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 2 of 2