

SLPI antibody

Cat. No. GTX39074

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
IHC-P	Assay dependent

Note : This antibody detects a band of approximately 14kDa in Daudi cell lysates.

Not tested in other applications.

Calculated MW 14 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 17 amino acid peptide located near the centre of human SLPI.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX39074 IHC-P Image

IHC-P analysis of mouse lung tissue using GTX39074 SLPI antibody. Heat mediated antigen retrieval using citrate buffer pH6.0.



GTX39074 WB Image

WB analysis of Daudi cell lysate using GTX39074 SLPI antibody at 1(A) and 2(B) ug/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).