

SNAP25 (cleaved) antibody [B371M]

Cat. No. GTX39119

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-Fr, ELISA, IHC
Reactivity	Human

References (3)

Package

500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
ELISA	1:100
IHC	Assay dependent

Not tested in other applications.

Product Note

This antibody recognizes the C. botulinum Toxin A cleaved SNAP-25 (183-197aa) in ELISA and cell staining assays. It does not cross-react with uncleaved SNAP-25, VAMP2, and Syntaxin 1.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide analogue SNAP-25 a.a. 183-197
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).