

## Salmonella species (A,B,C,D,E,F,G Groups) antibody [B345M]

## Cat. No. GTX39181

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	ICC/IF, ELISA, Sandwich ELISA
<b>Reactivity</b>	Salmonella

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX39181, Detection : GTX39180 / GTX39179**

**The detection antibody requires HRP conjugation. Please contact us for the conjugation service.**

Not tested in other applications.

## Product Note

Specific for common LPS core of Salmonella. Positive for all Salmonella O-serogroups tested: A, B, C1,C2, D, E1, E3, E4, F, G1 and G2. Negative for E. coli (mix), Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus and Legionella.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Mixture of 4 Salmonella serotypes (whole organisms).
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).