

Rhizopus arrhizus antibody [WSSA-RA-1]

Cat. No. GTX39191

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB, IHC-P, ELISA
Reactivity	Rhizopus arrhizus

References (1)

Package

250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1/50
ELISA	Assay dependent

Not tested in other applications.

Product Note Clone WSSA-RA-1 does not bind to water-soluble somatic antigens (WSSA) of Aspergillus spp.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Water-soluble somatic antigens (WSSA) from Rhizopus arrhizus.
Purification	Ammonium sulphate precipitation From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX39191 IHC-P Image

IHC analysis of renal tissue using GTX39191 *Rhizopus arrhizus* antibody [WSSA-RA-1].



For full product information, images and publications, please visit our [website](#).