

# Ricin, Ricinus communis agglutinin 60 (RCA60) B-chain antibody [RB999] (Biotin)

**Cat. No. GTX39237**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Application</b>	WB, ELISA, Sandwich ELISA
<b>Reactivity</b>	Ricinus communis

**Package**  
200 µg

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Sandwich ELISA : Capture: GTX39239, Detection: GTX39237**

Not tested in other applications.

**Product Note** Recognizes the B-chain of RCA60 and RCA120.

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	3.50 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Purification</b>	Protein G purified
<b>Conjugation</b>	Biotin

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).