

## Parathyroid Hormone antibody [A1/70], N-term

Cat. No. GTX39335

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC, RIA
Reactivity	Human

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC	2 µg/ml
RIA	20 ng/ml

Not tested in other applications.

**Product Note** Human PTH peptide (a.a. 15-25). Does not cross-react with Synthetic human PTH peptide (a.a. 1-3 or a.a. 1-10).

## Properties

Form	Liquid
Buffer	10mM Sodium Phosphate , 150mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic human PTH peptide (1-38) poly Lys. conjugated.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).