

Pancreatic alpha Amylase antibody

Cat. No. GTX39406

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Pig

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note

This antibody recognizes Amylase (Pancreatic), a digestive enzyme with a molecular weight of approximately 56 kDa secreted by the pancreas.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Porcine pancreatic alpha-amylase
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).