

## PSA (Total) antibody [906]

**Cat. No. GTX39433**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ELISA, Sandwich ELISA
<b>Reactivity</b>	Human

**Package**  
1 mg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX39441, Detection : GTX39433**

Not tested in other applications.

**Product Note** No cross reactivity with Human Albumin, AFP, CEA, PAP, CA-125, CA 19-9 or CA 15-3.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	15mM Potassium Phosphate, 0.85% NaCl
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Concentration</b>	7.50 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Human Total PSA
<b>Purification</b>	Purified by DEAE Chromatography
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).