

PSA (Free) antibody [981]

Cat. No. GTX39441

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX39441, Detection : GTX39433

Not tested in other applications.

Product Note No cross reactivity with Human Albumin, AFP, CEA, PAP, CA-125, CA 19-9 or CA 15-3.

Properties

Form	Liquid
Buffer	15mM Potassium Phosphate, 0.85% NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	6.41 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human PSA
Purification	Purified by DEAE Chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).