

Prolactin antibody [090-11710]

Cat No. GTX39483

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Application	ELISA, Sandwich ELISA
Reactivity	Human

Package
1 mg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

Calculated MW	26 kDa. (Note)
Product Note	Cross-reactivity: hGH <0.1%

PROPERTIES

Form	Liquid
Buffer	10mM Phosphate pH7.4, 150mM Sodium Chloride
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.57 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Source: Human Pituitary Gland
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).